



**Event COT102 Cotton**

**Real-time, Event-specific Polymerase Chain Reaction Method**

Syngenta Crop Protection, LLC  
3054 East Cornwallis Road  
PO Box 12257  
Research Triangle Park, NC 27709-2257, USA

For technical support contact [detection.method@syngenta.com](mailto:detection.method@syngenta.com)

## TABLE OF CONTENTS

|  |           |
|--|-----------|
| <b>TABLE OF CONTENTS</b>   | <b>2</b>  |
| <b>LIST OF TABLES</b>  | <b>3</b>  |
| <b>LIST OF FIGURES</b>   | <b>3</b>  |
| <b>LIST OF ACRONYMS AND ABBREVIATIONS</b>  | <b>4</b>  |
| <b>1.0 INTRODUCTION</b>  | <b>5</b>  |
| <b>2.0 METHOD</b>  | <b>5</b>  |
| 2.1 Principle of the COT102 Real-time, Event-specific Method.....                        | 5         |
| 2.2 Reagents and Equipment.....  | 6         |
| 2.3 DNA Concentration of Samples.....  | 7         |
| 2.4 Calibration Standards.....   | 7         |
| 2.5 The <i>sah7</i> -specific PCR Method.....  | 7         |
| 2.6 Event-specific PCR Method for the Detection and Quantification of<br>COT102 DNA..... | 9         |
| 2.7 Cycling Parameters.....  | 10        |
| 2.8 PCR Plate Layout.....  | 10        |
| 2.9 Data Analysis and Evaluation of Results.....   | 11        |
| <b>3.0 REFERENCES</b>  | <b>12</b> |

## LIST OF TABLES

|         |   |    |
|---------|---|----|
| TABLE 1 | Equipment and Materials .....   | 6  |
| TABLE 2 | Reagents, Buffers, and Solutions .....                                    | 6  |
| TABLE 3 | Dilution Scheme of the Calibration Standards .....                        | 7  |
| TABLE 4 | The <i>sah7</i> -specific PCR Primers and Probe Sequences .....           | 8  |
| TABLE 5 | Master Mix for the <i>sah7</i> -specific PCR Method .....                 | 9  |
| TABLE 6 | COT102 Real-time, Event-specific PCR Primers and Probe<br>Sequences ..... | 9  |
| TABLE 7 | Master Mix for the COT102-specific PCR Method .....                       | 10 |
| TABLE 8 | Cycling Parameters .....  | 10 |

## LIST OF FIGURES

|          |  |    |
|----------|--|----|
| FIGURE 1 | Sequence of the 115-bp Amplicon Generated by PCR<br>Amplification with the <i>sah7</i> -specific PCR Method from the A-<br>subgenome of Cotton ..... | 8  |
| FIGURE 2 | Sequence of the 123-bp Amplicon Generated by PCR<br>Amplification with the <i>sah7</i> -specific PCR Method from the D-<br>subgenome of Cotton ..... | 8  |
| FIGURE 3 | Sequence of the 101-bp Amplicon Generated by PCR<br>Amplification with the COT102 Real-time, Event-specific PCR<br>Method .....                      | 9  |
| FIGURE 4 | PCR Plate Layout .....   | 11 |

## LIST OF ACRONYMS AND ABBREVIATIONS

Definitions of International System of Units (SI) base units and derived units may be found in NIST (2011).

|       |  |
|-------|--|
| bp    | base pair  |
| CT    | cycle threshold  |
| DNA   | deoxyribonucleic acid  |
| DNase | deoxyribonuclease  |
| PCR   | polymerase chain reaction  |
| sah   | sinapis arabidopsis homolog  |
| TE    | tris (2-amino-2(hydroxymethyl)-propane-1,3-diol)/EDTA<br>(ethylenediaminetetraacetic acid) |

## 1.0 INTRODUCTION

Using the techniques of modern molecular biology, Syngenta has transformed cotton (*Gossypium hirsutum* L.) to produce Event COT102 cotton, a new cultivar that has controls against several lepidopteran pests of cotton.

A real-time, event-specific polymerase chain reaction (PCR) method was developed to detect and quantify COT102 deoxyribonucleic acid (DNA) extracted from seed and other plant material samples. The method consists of a cotton-specific PCR method as a reference, and an event-specific PCR method for detection of COT102 cotton.

## 2.0 METHOD

### 2.1 Principle of the COT102 Real-time, Event-specific Method

The COT102 real-time, event-specific PCR method is optimized for use on an Applied Biosystems<sup>®</sup> 7500 Fast Real-Time PCR System using the standard mode. The amount of PCR product is determined during each cycle (real-time) by measuring the fluorescence produced by a target-specific oligonucleotide probe labeled with two fluorescent dyes.

- COT102-specific PCR method: oligonucleotide probe with 6-FAM<sup>™</sup> as a reporter dye at its 5' end and TAMRA<sup>™</sup> as a quencher dye at its 3' end
- Cotton-specific PCR method: oligonucleotide probe with VIC<sup>®</sup> as a reporter dye at its 5' end and TAMRA<sup>™</sup> as a quencher dye at its 3' end

The 5'-nuclease activity of *Taq* DNA polymerase cleaves the probe and liberates the fluorescent moiety from proximity to the quenching moiety during the amplification process. The resulting increase in fluorescence during amplification is monitored and recorded.

A cotton-specific PCR method previously validated within the context of the validation of the event-specific quantitative method for Event 3006-210-23/281-23-235 cotton ([http://gmo-crl.jrc.ec.europa.eu/summaries/281-24-36\\_cotton\\_Protocol.pdf](http://gmo-crl.jrc.ec.europa.eu/summaries/281-24-36_cotton_Protocol.pdf), [http://gmo-crl.jrc.ec.europa.eu/summaries/3006-210-23\\_cotton\\_Protocol.pdf](http://gmo-crl.jrc.ec.europa.eu/summaries/3006-210-23_cotton_Protocol.pdf)) is used as a reference method.

For specific detection of COT102 genomic DNA, a 101-base pair (bp) fragment of the region that spans the 3' insert-to-plant genome junction in COT102 cotton is amplified using two specific primers. The forward primer binding site is located in the COT102 insert, the reverse primer binding site is located in the cotton genomic sequence, and the probe binding site spans the junction between the COT102 insert and the adjacent cotton genomic DNA.

## 2.2 Reagents and Equipment

All materials (e.g., vials, containers, pipette tips) should be suitable for PCR and molecular biology applications (Table 1). Materials should be deoxyribonuclease (DNase)-free, DNA-free, sterile, and unable to absorb protein or DNA. Table 2 contains a list of reagents, buffers, and solutions needed to perform the PCR method.

**TABLE 1 Equipment and Materials**

| Equipment and materials                            | Specification  |
|--|--|
| Applied Biosystems® 7500 Fast Real-Time PCR System | Applied Biosystems® part number 4351107 or equivalent                    |
| Vortex   | NeoLab VM-300 or equivalent  |
| Thermo-Fast® 96-well PCR plate                     | Abgene® catalog number AB-1400 or equivalent                             |
| Clear seal diamond heat sealing film               | Abgene® catalog number AB-0812 or equivalent                             |
| Thermo-sealer                                      | Abgene® catalog number AB-0384/240 or equivalent                         |
| Pipettes with adjustable volume                    | Eppendorf Research®, 2 – 20 µl, 20 – 200 µl, 100 – 1000 µl or equivalent |
| Filter tips  | Fitting the pipette models used  |
| Microcentrifuge tubes 1.5 ml                       | Roth catalog number 4182.1 or equivalent                                 |
| Microcentrifuge tubes 1.5 ml, screw lids           | Sarstedt part number 72.692.005 or equivalent                            |

**TABLE 2 Reagents, Buffers, and Solutions**

| Reagents, buffers, and solutions | Specification                                     |
|----------------------------------|---|
| TaqMan® Universal PCR Master Mix | Applied Biosystems® catalog number 4326708        |
| 1X TE buffer pH 8.0              | Applichem part number A2575, 1000                 |
| Water (HPLC Gradient Grade)      | Rotisolv® product number A511.1 or equivalent     |
| Salmon sperm DNA                 | Sigma-Aldrich® catalog number D7656 or equivalent |

TE = tris (2-amino-2(hydroxymethyl)-1,3-propanediol)/EDTA (ethylenediaminetetraacetic acid).  
HPLC = high performance liquid chromatography.

The handling of all reagents and controls should be carried out under sterile conditions and in a manner that precludes contamination of reagents or controls with exogenous DNA or undesired enzymatic activities (e.g., DNase).

The dilution buffer is 0.1X TE containing salmon sperm DNA. Resuspend 1 mg of salmon sperm DNA in 100 ml of 0.1 X TE (pH 8.0) to make a 10 ng DNA/ $\mu$ l 0.1X TE solution. Solution should be vortexed and stored at  $-20^{\circ}\text{C}$ .

### 2.3 DNA Concentration of Samples

For the unknown samples, the use of 100 ng of template DNA per reaction is recommended. This corresponds to approximately 42918 haploid copies of the cotton genome, assuming a genome weight of 2.33 pg (Arumuganathan and Earle 1991).

### 2.4 Calibration Standards

The method format uses five calibration standards containing different amounts of COT102 genomic DNA (Standard 1 through Standard 5); each calibration standard is analyzed in both the COT102-specific and *Sinapis Arabidopsis* Homolog (*sah*) 7-specific PCR methods. The calibration standards are produced by preparing solutions of 20 ng/ $\mu$ l (100 ng/reaction) of total genomic DNA with 10%, 5%, 1%, 0.5%, and 0.07% copy/copy COT102 DNA mixed with nontransgenic cotton genomic DNA. Table 3 lists the total DNA content, as well as the total copy numbers of COT102, for each calibration standard. COT102 copy numbers are based on the fact that the COT102 seeds are homozygous.

**TABLE 3 Dilution Scheme of the Calibration Standards**

| Sample            | Standard 1 | Standard 2 | Standard 3 | Standard 4 | Standard 5 |
|-------------------|------------|------------|------------|------------|------------|
| Total DNA content | 100 ng     | 100 ng     | 100 ng     | 100 ng     | 100 ng     |
| COT102 DNA copies | 4292       | 2146       | 429        | 215        | 30         |

### 2.5 The *sah7*-specific PCR Method

For relative quantification of COT102 DNA, two specific primers (Table 4) are used to amplify a fragment of the cotton endogenous gene *sah7*. This cotton-specific PCR method has been validated by the JRC ([http://gmo-crl.jrc.ec.europa.eu/summaries/281-24-36\\_cotton\\_Protocol.pdf](http://gmo-crl.jrc.ec.europa.eu/summaries/281-24-36_cotton_Protocol.pdf), [http://gmo-crl.jrc.ec.europa.eu/summaries/3006-210-23\\_cotton\\_Protocol.pdf](http://gmo-crl.jrc.ec.europa.eu/summaries/3006-210-23_cotton_Protocol.pdf)). Whereas primer and probe sequences and concentrations were taken from the original protocol, the buffer system was modified and validated. The gene *sah7* is present in both the A- and the D-subgenomes of *G. hirsutum*. The A-subgenome specific copy differs from the D-subgenome specific copy by several single or double nucleotide substitutions and small deletions or insertions. Primers and probe sequences of the cotton-specific reference PCR system match perfectly without any single mismatch to both subgenomes' gene copies. However, due to sequence deviations within the amplified region, the sizes of the amplicons resulting from the A- and D- subgenomes differ slightly, being 115 bp and 123 bp in length (Figures 1 and 2), respectively. The amount of PCR product is determined during each cycle in real-time by measuring

fluorescence produced by a *sah7*-specific oligonucleotide probe (Table 4) labeled with VIC<sup>®</sup> as a reporter dye at its 5' end and TAMRA<sup>™</sup> as a quencher at its 3' end.

**TABLE 4 The *sah7*-specific PCR Primers and Probe Sequences**

| Primer/probe name  | Primer sequence 5' to 3'   |
|--------------------|--|
| SAH7-uni-f1 primer | AGTTTGTAGGTTTTGATGTTACATTGAG   |
| SAH7-uni-r1 primer | GCACTCTTGAACCGCCTACTG  |
| SAH7-uni-s1 probe  | VIC <sup>®</sup> -AAACATAAAATAATGGGAACAACCATGACATGT-TAMRA <sup>™</sup> |

**FIGURE 1 Sequence of the 115-bp Amplicon Generated by PCR Amplification with the *sah7*-specific PCR Method from the A-subgenome of Cotton**

SAH7-uni-f1 primer SAH7-uni-s1 probe  
AGTTTGTAGGTTTTGATGTTACATTGAGTGACAGTGAATGAAAGGGTGTGTAAACATAAAATAATGGGAA  
CAACCATGACATGTTGGACTGGATCAGTAGGCGGTTCAAAGATGC  
SAH7-uni-r1 primer

Primer and probe binding sites are underlined.

**FIGURE 2 Sequence of the 123-bp Amplicon Generated by PCR Amplification with the *sah7*-specific PCR Method from the D-subgenome of Cotton**

SAH7-uni-f1 primer SAH7-uni-s1 probe  
AGTTTGTAGGTTTTGATGTTACATTGAGTGAAAGCGAATGATAGGGTGTCTAAACATAAAATAATGGGAA  
CAACCATGACATGTTGGACTGGAGGTCAAATCAGTAGGCGGTTCAAAGATGC  
SAH7-uni-r1 primer

Primer and probe binding sites are underlined.

All reagents should be thawed, as necessary, and thoroughly mixed before each use. The master mix for the *sah7*-specific PCR method that contains all components of the PCR method except template DNA (Table 5) can be prepared in sufficient quantities before reactions are performed.





**TABLE 7 Master Mix for the COT102-specific PCR Method**

| Components                       | Initial concentration | Final concentration | µl/reaction  |
|----------------------------------|-----------------------|---------------------|--------------|
| TaqMan® Universal PCR Master Mix | 2X                    | 1X                  | 12.50        |
| COT102_3_89F primer              | 10 µM                 | 600 nM              | 1.50         |
| COT102_3_181R primer             | 10 µM                 | 600 nM              | 1.50         |
| COT102_3_115T probe              | 10 µM                 | 150 nM              | 0.375        |
| Nuclease-free water              | N/A                   | N/A                 | 4.125        |
| <b>Total volume<sup>a</sup></b>  | N/A                   | N/A                 | <b>20.00</b> |

<sup>a</sup>Total PCR reaction is 25 µl (20 µl master mix and 5 µl genomic DNA template).  
N/A = not applicable.

## 2.7 Cycling Parameters

An Applied Biosystems® 7500 Fast Real-Time PCR System using the standard mode in the instrument set-up is used with the COT102 real-time, event-specific PCR method and the *sah7*-specific PCR method. Table 8 shows the PCR cycling conditions.

**TABLE 8 Cycling Parameters**

| Cycle | Step | Temperature (°C) | Time (seconds) | Data acquisition | Number of cycles |
|-------|------|------------------|----------------|------------------|------------------|
| A     | 1    | 50               | 120            | No               | 1                |
| B     | 1    | 95               | 600            | No               | 1                |
| C     | 1    | 95               | 15             | No               | 45               |
|       | 2    | 60               | 60             | Yes              |                  |

## 2.8 PCR Plate Layout

The calibration of the COT102 real-time, event-specific method includes five calibration standards, each containing different percentages of COT102 cotton genomic DNA, with a total amount of 100 ng cotton DNA. The COT102 DNA content of the standard samples ranges from 0.07% to 10% copy/copy. Each calibration standard is analyzed in duplicate with the COT102-specific PCR method, as well as with the *sah7*-specific PCR method. Duplicate negative control samples (1X TE buffer) containing no template DNA verify the purity of the reagents. Unknown samples are analyzed in triplicate for both the COT102-specific and *sah7*-specific PCR methods. Figure 4 summarizes the PCR plate layout.

**FIGURE 4 PCR Plate Layout**

|   | 1          | 2 | 3          | 4          | 5          | 6 | 7          | 8 | 9                | 10         | 11               | 12 |
|---|------------|---|------------|------------|------------|---|------------|---|------------------|------------|------------------|----|
| A | Standard 1 |   | Standard 1 |            | Standard 2 |   | Standard 2 |   | Standard 3       |            | Standard 3       |    |
| B | Standard 4 |   | Standard 4 |            | Standard 5 |   | Standard 5 |   | Negative control |            | Negative control |    |
| C | Unknown 1  |   |            | Unknown 1  |            |   | Unknown 2  |   |                  | Unknown 2  |                  |    |
| D | Unknown 3  |   |            | Unknown 3  |            |   | Unknown 4  |   |                  | Unknown 4  |                  |    |
| E | Unknown 5  |   |            | Unknown 5  |            |   | Unknown 6  |   |                  | Unknown 6  |                  |    |
| F | Unknown 7  |   |            | Unknown 7  |            |   | Unknown 8  |   |                  | Unknown 8  |                  |    |
| G | Unknown 9  |   |            | Unknown 9  |            |   | Unknown 10 |   |                  | Unknown 10 |                  |    |
| H | Unknown 11 |   |            | Unknown 11 |            |   | Unknown 12 |   |                  | Unknown 12 |                  |    |

**Standard 1 to 5:** calibration standards.

**Unknown 1 to 12:** unknown samples.

**Negative control:** samples containing no template DNA (1X TE buffer).

**White:** *sah7*-specific PCR method.

**Gray:** COT102-specific PCR method.

## 2.9 Data Analysis and Evaluation of Results

Analyze the results according to the following procedure:

- a) Place the threshold in the region of exponential amplification across all of the amplification plots. This region is recognized in the logarithmic view of the amplification plots as the portion of the plot that is linear. Do not place the threshold line in the plateau phase or in the initial linear phase of amplification. Place the threshold line clearly above the background fluorescence and above the level where splitting or fork effects between replicates can be observed.
- b) Determine the cycle number at which the threshold line crosses the first amplification curve and set the baseline three cycles before that value (e.g., if the earliest cycle threshold [CT] = 24, set the baseline crossing at 21, or CT of 24 – 3).
- c) Save the settings.
- d) Repeat the procedure described in a) and b) on the amplification plots of the other method.
- e) Save the settings and export all the data into a Microsoft® Office Excel® file for further calculations.

A calibration curve can be produced by plotting mean delta CT ( $\Delta$ CT) values of calibration samples against the logarithm of the respective percent COT102 DNA contents; the slope (a) and the intercept (b) of the calibration curve ( $y = ax + b$ ) can then be used to calculate the percent COT102 DNA content of the reference samples based on their normalized mean  $\Delta$ CT values.

### 3.0 REFERENCES

Arumuganathan K, Earle ED. 1991. Nuclear DNA content of some important plant species. *Plant Mol Biol Rep* 9:208–218.

NIST. 2011. *The NIST Reference on Constants, Units, and Uncertainty*. Gaithersburg, MD: National Institute of Standards and Technology. <http://physics.nist.gov/cuu> (updated June 2, 2011).